

DC Medical Equipment, LLC

Safety Precautions/Contraindications For Use:

- The stim pads given have a lifespan of 3-5 uses, we recommend switching them out after 5 uses and disposing of the used pads
- After use, return electrodes to original packaging to prolong usage and safety
- Do not place electrodes over cuts or wounds (or areas with a skin condition such as Eczema or Psoriasis)
- Do not place electrodes over the heart or directly over the spinal cord
- Do not place electrodes on front of neck, face or head
- Do not use if you have a pacemaker
- Do not use if you are pregnant
- Only use the unit for the specific use prescribed by your physician
- Do not allow anyone else to use or borrow your unit